

SERVOX® DIGITAL ELECTROLARYNX

REF 70000



- ▶ A modern, high-quality electric speaking aid
- ▶ The speaking aid is applied to the user's throat and generates sound oscillations, which are transmitted to the pharyngeal, oral and nasal area through a membrane.
- ▶ These oscillations are converted through natural speaking motions and form an understandable replacement voice.
- ▶ The manual basic settings allow customisation and saving of the audio frequency and sound qualities. Optional software is available for further customisation options.
- ▶ The volume control can be adjusted individually to suit the user's needs using a rotary knob.
- ▶ A soft screw cap for a more muted sound is available optionally.

TECHNICAL DATA

Volume	max. 87 dB
Dimensions	L 128 x Ø 35 mm
Weight	approx. 110 g excl. batteries
Talk time*	4-5 hours*
Voltage	1.2 V/1.5 V
Battery type:	AAA NiMH (HR03)

*The talk time depends on the selected battery capacity. The battery life specified by the manufacturer must be observed.

SCOPE OF DELIVERY

SERVOX® digital Electrolarynx REF 70000 Sprechhilfe

1 x SERVVOX® digital Electrolarynx speaking aid
1 x SERVVOX® digital Electrolarynx screw cap hard tone
1x SERVVOX® digital Electrolarynx strap set
1x SERVVOX® digital Electrolarynx belt pouch
1x SERVVOX® digital Electrolarynx battery holder
4x Rechargeable battery (type AAA)
1x Charger Varta LCD Plug
1x Instructions for use

Optionally available:

SERVOX® digital Electrolarynx REF 70001 Screw cap hard tone

SERVOX® digital Electrolarynx REF 70002 Screw cap soft tone

SERVOX® digital Electrolarynx REF 70003 Strap set

SERVOX® digital Electrolarynx REF 70004 Oral tube with adapter

1x Adapter
2x Oral tube short
2x Oral tube long

SERVOX® digital Electrolarynx REF 70005 Battery holder

Travel adapter

CE 0482